

Slit Knife

Disposable ophthalmic knives

NO.1

Slit Knife



Introduction


Blade width from 1.0mm to 3.5mm, with or without protective shield.

Straight Knife

Disposable ophthalmic knives

NO.2

Straight Knife



Introduction

Blade width 1.6mm with 5 different angles, with or without protective shield.

NO.	without protection shield	with protection shield	Angle		Blade Width
Slit Knife			$\alpha \pm 5^\circ$	$\alpha \pm 5^\circ$	$w \pm 0.05(\text{mm})$
1	E3SLU1-1.0	E3SLU2-1.0	43°	45°	1.0
2	E3SLU1-1.1	E3SLU2-1.1			1.1
3	E3SLU1-1.2	E3SLU2-1.2			1.2
4	E3SLU1-1.3	E3SLU2-1.3			1.3
5	E3SLU1-1.4	E3SLU2-1.4			1.4
6	E3SLU1-1.6	E3SLU2-1.6			1.6
7	E3SLU1-1.8	E3SLU2-1.8			1.8
8	E3SLU1-1.9	E3SLU2-1.9	45°	45°	8
9	E3SLU1-2.0	E3SLU2-2.0			9
10	E3SLU1-2.1	E3SLU2-2.1			10
11	E3SLU1-2.2	E3SLU2-2.2			11
12	E3SLU1-2.3	E3SLU2-2.3			12
13	E3SLU1-2.4	E3SLU2-2.4			13
14	E3SLU1-2.5	E3SLU2-2.5			14
15	E3SLU1-2.6	E3SLU2-2.6			15
16	E3SLU1-2.65	E3SLU2-2.65			16
17	E3SLU1-2.7	E3SLU2-2.7			17
18	E3SLU1-2.75	E3SLU1-2.75	47°	45°	18
19	E3SLU1-2.8	E3SLU1-2.8			19
20	E3SLU1-2.85	E3SLU1-2.85			20
21	E3SLU1-2.9	E3SLU1-2.9			21
22	E3SLU1-3.0	E3SLU1-3.0			22
23	E3SLU1-3.2	E3SLU1-3.2			23
24	E3SLU1-3.5	E3SLU1-3.5			24
25	E3SLS1-1.0	E3SLS1-1.0	43°	Straight	25
26	E3SLS1-2.8	E3SLS1-2.8			26
27	E3SLS1-3.0	E3SLS1-3.0	47°	Straight	27
Straight Knife			$\alpha \pm 3^\circ$	αI	$w \pm 0.05 (\text{mm})$
28	E3STS1-0°	E3STS2-0°	0°	Straight	1.6
29	E3STS1-15°	E3STS2-15°	15°		
30	E3STS1-22.5°	E3STS2-22.5°	22.5°		
31	E3STS1-30°	E3STS2-30°	30°		
32	E3STS1-45°	E3STS2-45°	45°		