


Crescent Knife

Disposable ophthalmic knives

NO.3

Crescent Knife



Introduction


Blade width from 2.0mm to 2.8mm, with or without protective shield.

MVR Knife

Disposable ophthalmic knives

NO.4

MVR Knife



Introduction

Blade width from 0.7mm to 15mm, with or without protective shield.

Crescent Knife			αl ±5°		w±0.1 (mm)
33	E3CRU1-2.0	E3CRU2-2.0	60°		2.0
34	E3CRU1-2.3	E3CRU2-2.3			2.3
35	E3CRU1-2.5	E3CRU2-2.5			2.5
36	E3CRU1-2.6	E3CRU2-2.6			2.6
37	E3CRU1-2.8	E3CRU2-2.8			2.8
38	E3CRD1-2.3	E3CRD2-2.3			2.3
39	E3CRS1-2.0	E3CRS2-2.0	Straight		2.0
40	E3CRS1-2.3	E3CRS2-2.3	Straight		2.3
MVR Knife – I			α±5 °	αl±5 °	w±0.05 (mm)
41	E3MRU1-19g	E3MRU2-19g	25°	45°	1.5
42	E3MRU1-20g	E3MRU2-20g	25°		1.2
43	E3MRU1-21g	E3MRU2-21g	25°		1.1
44	E3MRU1-23g	E3MRU2-23g	32°		0.9
45	E3MRU1-24g	E3MRU2-24g	32°		0.8
46	E3MRU1-23.5g	E3MRU2-23.5g	32°		0.7
47	E3MRS1-19g	-	25°		1.5
48	E3MRS1-20g	-	25°		1.2
49	E3MRS1-23g	-	32°		0.9
50	E3MRS1-24g	-	32°		0.8
R Slit Knife– I			α	αl±5 °	w±0.1 (mm)
51	E3RSU1-2.6	E3RSU2-2.6	--	45 °	2.6
52	E3RSU1-3.2	E3RSU2-3.2			3.2
53	E3RSU1-3.5	E3RSU2-3.5			3.5
54	E3RSU1-3.7	E3RSU2-3.7			3.7
55	E3RSU1-4.1	E3RSU2-4.1			4.1
56	E3RSU1-5.2	E3RSU2-5.2			5.2
57	E3RSU1-5.5	E3RSU2-5.5			5.5
Pterygium Knife – I			α±3 °	αl	w±0.1 (mm)
58	E3PMS1-2.4	E3PMS2-2.4	--	Straight	2.4