


R Slit Knife

Disposable ophthalmic knives

NO.5

R Slit Knife



Introduction

Blade width from 2.5mm to 5.5mm, with or without protective shield.

Pterygium Knife

Disposable ophthalmic knives

NO.6

Pterygium Knife



Introduction

Blade width 2.4mm, with or without protective shield.

Crescent Knife			$\alpha \pm 5^\circ$		w ± 0.1 (mm)
33	E3CRU1-2.0	E3CRU2-2.0	60°		2.0
34	E3CRU1-2.3	E3CRU2-2.3			2.3
35	E3CRU1-2.5	E3CRU2-2.5			2.5
36	E3CRU1-2.6	E3CRU2-2.6			2.6
37	E3CRU1-2.8	E3CRU2-2.8			2.8
38	E3CRD1-2.3	E3CRD2-2.3			2.3
39	E3CRS1-2.0	E3CRS2-2.0	Straight		2.0
40	E3CRS1-2.3	E3CRS2-2.3	Straight		2.3
MVR Knife - I			$\alpha \pm 5^\circ$	$\alpha \pm 5^\circ$	w ± 0.05 (mm)
41	E3MRU1-19g	E3MRU2-19g	25°	45°	1.5
42	E3MRU1-20g	E3MRU2-20g	25°		1.2
43	E3MRU1-21g	E3MRU2-21g	25°		1.1
44	E3MRU1-23g	E3MRU2-23g	32°		0.9
45	E3MRU1-24g	E3MRU2-24g	32°		0.8
46	E3MRU1-23.5g	E3MRU2-23.5g	32°		0.7
47	E3MRS1-19g	-	25°		1.5
48	E3MRS1-20g	-	25°		1.2
49	E3MRS1-23g	-	32°		0.9
50	E3MRS1-24g	-	32°		0.8
R Slit Knife- I			α	$\alpha \pm 5^\circ$	w ± 0.1 (mm)
51	E3RSU1-2.6	E3RSU2-2.6	--	45 °	2.6
52	E3RSU1-3.2	E3RSU2-3.2			3.2
53	E3RSU1-3.5	E3RSU2-3.5			3.5
54	E3RSU1-3.7	E3RSU2-3.7			3.7
55	E3RSU1-4.1	E3RSU2-4.1			4.1
56	E3RSU1-5.2	E3RSU2-5.2			5.2
57	E3RSU1-5.5	E3RSU2-5.5			5.5
Pterygium Knife - I			$\alpha \pm 3^\circ$	α	w ± 0.1 (mm)
58	E3PMS1-2.4	E3PMS2-2.4	--	Straight	2.4