



TRIVA

TRIVA-αAY YELLOW
TRIVAT-αAY YELLOW



TRIVA-αAY YELLOW	SL	TRIVAT-αAY YELLOW	SL
--	----	---	----



Type	Trifocal posterior chamber IOL, one-piece, foldable blue-light protection	Trifocal-toric posterior chamber IOL, one-piece, foldable blue-light protection
Optic diameter	6.0 mm	6.0 mm
Total diameter	12.5 mm	12.5 mm
Material	Hydrophilic acrylic with UV-absorber Glistening-free	Hydrophilic acrylic with UV-absorber Glistening-free
Optic features	Central diffractive aspheric anterior surface with a refractive optic periphery Aberration-free 360° lens epithelial cell barrier Intermediate addition +1.75 D Near addition +3.5 D	Central diffractive aspheric anterior surface with a refractive optic periphery Aberration-free 360° lens epithelial cell barrier Intermediate addition +1.75 D Near addition +3.5 D
Haptic design	C-loop	C-loop
A-constants	118.1 (ultrasonic measurement) 118.4 (optical measurement)	
Diopter range	+10.0 to +30.0 in 0.5 D steps	SE: +10.0 to +30.0 in 0.5 D steps Cyl.: 1.0 to 6.0 in 0.5 D steps

Surgical

Type	Product Model	Diopter Range	Optical Diameter	Haptic Design	Aspherical	Toric	Blue light protection	Preloaded
Trifocal	TRIVA-aAY YELLOW	+10.0~+30.0D in 0.5D step	6.0 mm	C-loop	●		●	Optional
Trifocal	TRIVAT-aAY YELLOW	SE: +10.0~+30.0D in 0.5D step Cyl.: 1.0~6.0D in 0.5D step	6.0 mm	C-loop	●	●	●	Optional
Bifocal	DiffRACTIVA-aAY YELLOW	+10.0D to +30.0D (Biconvex)	6.0 mm	C-loop	●		●	●
Bifocal	DiffRACTIVA-aA							
Monofocal	TORICA-aAY YELLOW	SE: -20.0~+60.0D in 0.5D step Cyl.: 1.0~20.0D in 0.5D step	6.0 mm	C-loop	●	●	●	Optional
Monofocal	TORICA-aA							
Monofocal	ASPIRA-aXAY YELLOW	Preloaded in SAFELOADER® +10.0~+30.0D in 0.5D step	7.0 mm	Cut-out haptics	●		●	Optional
Monofocal	ASPIRA-aXA	In COMPACT LINE -10.0~+30.0D in 1.0D step +10.0~+30.0D in 0.5D step						
Monofocal	ASPIRA-aAY YELLOW	Preloaded in SAFELOADER® 10.0 bis 30.0 in 0.5 D steps	6.0 mm	C-loop	●		●	Optional
Monofocal	ASPIRA-aA	In COMPACT LINE -20,0 to 60.0 in 1.0 D steps 10.0 to 30.0 in 0.5 D steps						
Monofocal	MC XII ASP	-6.0D~+40.0 D	5.5-7.0 mm	H-beam	●		●	