


Diffractiva®

Diffractiva®-aAY YELLOW
Diffractiva®-aA



Surgical

	DIFFRACTIVA®-aAY YELLOW SL	DIFFRACTIVA®-aA SL
--	--	---------------------------------



Type	Multifocal IOL Bluelight filtering	Multifocal IOL
Optic diameter	6.0 mm	7.0mm, 6.5mm
Total diameter	12.5 mm	11.0mm
Material	Hydrophilic, Glistening-free acrylic, UV absorber	
Optic features	Aspheric anterior surface with central diffractive zone +3.5D add power 360° LEC barrier on posterior surface One piece design with C loop	
Spherical Aberration	Aberration-free	Aberration-free
Diopter range	+10.0 D ~ +30.0 D in 0.5D	
Injector system	SAFELOADER®	SAFELOADER®
A-constants	118.3 (ultrasonic measurement) 118.6 (optical measurement)	
XS Incision	1.8mm	1.8mm

Surgical

Type	Product Model	Diopter Range	Optical Diameter	Haptic Design	Aspherical	Toric	Blue light protection	Preloaded
Trifocal	TRIVA-aAY YELLOW	+10.0~+30.0D in 0.5D step	6.0 mm	C-loop	●		●	Optional
Trifocal	TRIVAT-aAY YELLOW	SE: +10.0~+30.0D in 0.5D step Cyl.: 1.0~6.0D in 0.5D step	6.0 mm	C-loop	●	●	●	Optional
Bifocal	DiffRACTIVA-aAY YELLOW	+10.0D to +30.0D (Biconvex)	6.0 mm	C-loop	●		●	●
Bifocal	DiffRACTIVA-aA							
Monofocal	TORICA-aAY YELLOW	SE: -20.0~+60.0D in 0.5D step Cyl.: 1.0~20.0D in 0.5D step	6.0 mm	C-loop	●	●	●	Optional
Monofocal	TORICA-aA							
Monofocal	ASPIRA-aXAY YELLOW	Preloaded in SAFELOADER® +10.0~+30.0D in 0.5D step	7.0 mm	Cut-out haptics	●		●	Optional
Monofocal	ASPIRA-aXA	In COMPACT LINE -10.0~+30.0D in 1.0D step +10.0~+30.0D in 0.5D step						
Monofocal	ASPIRA-aAY YELLOW	Preloaded in SAFELOADER® 10.0 bis 30.0 in 0.5 D steps	6.0 mm	C-loop	●		●	Optional
Monofocal	ASPIRA-aA	In COMPACT LINE -20,0 to 60.0 in 1.0 D steps 10.0 to 30.0 in 0.5 D steps						
Monofocal	MC XII ASP	-6.0D~+40.0 D	5.5-7.0 mm	H-beam	●		●	