



Aspira®

Aspira®-αXAY Yellow
Aspira®-αXA



Surgical

ASPIRA®-αXAY YELLOW SL	ASPIRA®-αXA SL
--	--



Type	Monofocal posterior chamber IOL, one-piece, foldable, blue-light protection optional	Monofocal posterior chamber IOL, one-piece, foldable, blue-light protection optional
Optic diameter	7.0 mm	7.0 mm
Total diameter	11.0 mm	11.0 mm
Material	Hydrophilic acrylic with UV-absorber, glistening-free	
Optic features	XL optic, aspheric anterior surface, aberration-free, 360° lens epithelial cell barrier	XL optic, aspheric anterior surface, aberration-free, 360° lens epithelial cell barrier
Haptic design	Cut-out haptics	Cut-out haptics
A-constants	118.1 (ultrasonic measurement) 118.4 (optical measurement)	

Diopter range	Preloaded in SAFELOADER® +10.0~+30.0D in 0.5D step	Preloaded in SAFELOADER® +10.0~+30.0D in 0.5D step
	In COMPACT LINE -10.0~+30.0D in 1.0D step +10.0~+30.0D in 0.5D step	In COMPACT LINE -10.0~+30.0D in 1.0D step +10.0~+30.0D in 0.5D step

Surgical

Type	Product Model	Diopter Range	Optical Diameter	Haptic Design	Aspherical	Toric	Blue light protection	Preloaded
Trifocal	TRIVA-aAY YELLOW	+10.0~+30.0D in 0.5D step	6.0 mm	C-loop	●		●	Optional
Trifocal	TRIVAT-aAY YELLOW	SE: +10.0~+30.0D in 0.5D step Cyl.: 1.0~6.0D in 0.5D step	6.0 mm	C-loop	●	●	●	Optional
Bifocal	Diffraction-aAY YELLOW	+10.0D to +30.0D (Biconvex)	6.0 mm	C-loop	●		●	●
Bifocal	Diffraction-aA							
Monofocal	TORICA-aAY YELLOW	SE: -20.0~+60.0D in 0.5D step Cyl.: 1.0~20.0D in 0.5D step	6.0 mm	C-loop	●	●	●	Optional
Monofocal	TORICA-aA							
Monofocal	ASPIRA-aXAY YELLOW	Preloaded in SAFELOADER® +10.0~+30.0D in 0.5D step	7.0 mm	Cut-out haptics	●		●	Optional
Monofocal	ASPIRA-aXA	In COMPACT LINE -10.0~+30.0D in 1.0D step +10.0~+30.0D in 0.5D step						
Monofocal	ASPIRA-aAY YELLOW	Preloaded in SAFELOADER® 10.0 bis 30.0 in 0.5 D steps	6.0 mm	C-loop	●		●	Optional
Monofocal	ASPIRA-aA	In COMPACT LINE -20.0 to 60.0 in 1.0 D steps 10.0 to 30.0 in 0.5 D steps						
Monofocal	MC XII ASP	-6.0D~+40.0 D	5.5-7.0 mm	H-beam	●		●	